MD

A detailed line drawing of the L3000 system components. The diagram includes a coiled cable (1) with a connector (2) at one end and a control unit (3) in the middle. A long, thin probe (4) is attached to the control unit. A separate component (5) is shown with a button labeled 'L3000'. Other callouts point to specific parts: 6 points to the 'L3000' button, 7 points to the control unit, 8 points to the cable, and 9 points to the coiled cable.

A close-up photograph of the cryoICE BOX interface. The device is a grey rectangular box with the text "cryoICE BOX" printed on the top right. On the left side, there is a heart icon and a small symbol. Three numbered callouts point to specific components: callout 2 points to a circular port with a cable; callout 7 points to a cable connector; and callout 8 points to a metal fitting on a cable. The cables are connected to the device, and the overall setup is part of a medical or research equipment.

en 13 mm
(0.5 inch)
vi 13 mm
(0,5 inch)



en 26 mm
(1 inch)
vi 26 mm
(1 inch)

en 13 mm
(0.5 inch)
vi 13 mm
(0,5 inch)

A diagram showing a curved pipe section. The pipe is represented by two parallel lines forming an arc. A vertical dimension line on the right indicates a height of 5 cm (2.0 inches) from the horizontal baseline to the top of the pipe section.

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en

Operational	Storage	Transit
Temperature: 10°C/50°F to 40°C/104°F	Temperature: -29°C/-20°F to 60°C/140°F	Temperature: -29°C/-20°F to 60°C/140°F
Relative Humidity: 15% to 90%	Relative Humidity: 30% to 85%	Relative Humidity: 30% to 85%
Atmospheric Pressure: 98 to 105 kPa (14.2 to 15.2 psi)	Atmospheric Pressure: 98 to 105 kPa (14.2 to 15.2 psi)	N/A

[5] Malleable Section of PROBE

The ACM components are not suitable for use in the presence of a flammable anesthetic mixture which can cause a fire or explosion, resulting in user and patient injury or death.

If the sterile package is dropped and/or damaged or the sterile barrier is breached, discard device and DO NOT USE.
Breach of sterile barrier can lead to infection.

NOTE: If the same bend is desired in a different plane, do not twist the Malleable Section of PROBE; re-straighten the Malleable Section of PROBE and create the same bend in the desired plane.

Do not bend Malleable Section of the PROBE during FREEZE or DEFROST mode. It can cause a high pressurized gas leak that can potentially lead to tissue perforation, unintended damage, or injury to user.

⚠ CAUTION: The Malleable Section of PROBE has a limited functional life; if greater than 8 bends are intended, it is recommended to use a second probe.

- AtriCure®**
6. It is always recommended to use the TOOL to create desired bends. The TOOL has two ends, the smaller end radius is 13 mm (0.5 inches) and the larger end radius is 26 mm (1.0 inch).
7. Typical procedures may require the following bend profiles created with the use of the TOOL, as illustrated in **Figure 4**.

FORMING THE RIGID PROBE SHAFT

⚠ CAUTION: The distal end of the Rigid PROBE Shaft should not be bent more than 5 cm (2.0 inches) from straight, as illustrated in **Figure 5**.

⚠ CAUTION: Do not use the PROBE if damaged as it may result in device malfunction. Repetitive bends in the same location could cause damage to the Rigid PROBE Shaft. The PROBE has a limited functional life; if greater than 7 Rigid Probe Shaft bend cycles are intended, it is recommended to use a second probe.

USING THE PROBE TO PERFORM CRYOABLATION

NOTE: The PROBE ablates tissue via cryogenic energy delivered to the Malleable Section of PROBE. Cryoadhesion of the Malleable Section of PROBE to tissue can occur when the PROBE reaches a temperature of 0°C (32°F) or below. Other portions of the PROBE, including the Rigid PROBE Shaft, can become cold, and should be handled with appropriate care.

⚠ CAUTION: Do not use the PROBE if damaged as it may result in device malfunction. The PROBE has a limited functional life; if greater than 14 Freeze/Defrost cycles are intended, it is recommended to use a second probe.

- With the PROBE in air, prime the system with a Pre-Freeze cycle: Set the CONSOLE Ablation Timer to 30 seconds and press the Activation Button to engage FREEZE Mode. Wait for the system to cycle through FREEZE and DEFROST, or manually advance via the Activation Button.
- During the FREEZE cycle, if there are leaks in the blue inlet/orange exhaust connector, the sound of gas leaking will be heard and/or frost will appear on the connections. Replace the device before continuing with the procedure.

 WARNING

- Set the Ablation Timer to the desired ablation time. The timer is preset to a default of 120 seconds.
- Navigate the PROBE to the target ablation site:
 - Identify the target site.
 - Reach the Malleable Section of PROBE through an appropriate-sized incision to the target.
 - Under direct visualization, place the Malleable Section of PROBE against the target tissue.

 WARNING

- Using the Retractable Handle, apply gentle pressure to the Malleable Section of PROBE, and avoid any PROBE movement until after the freeze cycle completes.
- Under direct visualization ensure that the PROBE Malleable Section of PROBE and Rigid PROBE Shaft are not in contact with other anatomical structures not intended for ablation.

 WARNING

- Press the activation button or use the optional ACM footswitch to engage FREEZE Mode for the desired length of time. The system will automatically cycle from FREEZE to DEFROST after the Ablation Timer has expired.

 WARNING

- Wait until the PROBE temperature has warmed to above 0°C (32°F) before attempting to remove the Malleable Section of PROBE from the ablation site or moving the Rigid PROBE Shaft.
- ⚠ CAUTION:** Use care while the CONSOLE is in Defrost Mode, as during N₂O gas venting, the PROBE may cool sufficiently to cause cryoadhesion.

NOTE: If PROBE does not reach desired DEFROST temperature, apply warm, sterile, saline to the tissue and PROBE area as necessary.

- After the CONSOLE is in Ready Mode and the PROBE temperature is above 0°C (32°F), repeat steps (11) to (16) to create additional cryoablation lesions.

DISCONNECTING AND DISPOSING OF THE PROBE

- Close N₂O Cylinder by turning the Valve fully clockwise.
- Pull the red N₂O Manual Exhaust Knob or press the N₂O Exhaust Switch on the back of the CONSOLE to fully depressurize the system.

⚠ CAUTION: Ensure the CONSOLE is in Ready Mode before attempting to disconnect the PROBE. The sudden release of pressurized gas may cause the PROBE to recoil, which may injure the operator or patient.

- Disconnect the PROBE from the CONSOLE and discard.

 WARNING

DISPOSAL

After use this device should be treated as medical waste and disposed of following hospital protocol.

SERIOUS INCIDENT

Any serious incident that has occurred in relation to this device should be reported to AtriCure and the competent authority of the Member State in which the user and/or patient is located.

RETURN OF USED PRODUCT

If for any reason this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping. If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

TROUBLESHOOTING

PROBLEM	POTENTIAL CAUSE	SOLUTION
PROBE does not reach desired defrost temperature after freeze.	Plugged gas supply path.	Manually defrost by applying warm saline to tissue and probe as necessary.
PROBE does not reach the proper temperature.	Empty or low N ₂ O cylinder.	Replace low or empty N ₂ O cylinder.
	Gas not flowing, tubing is restricted.	Verify PROBE tubing is not pinched.
	Gas leak in PROBE Shaft or tubing.	Replace PROBE.
	N ₂ O tank valve closed.	Fully open N ₂ O tank valve.
CONSOLE displays “----”.	Thermocouple Connectors not fully plugged into the CONSOLE.	Plug Thermocouple Connectors all the way into the CONSOLE ports.
	PROBE internal wires are broken.	Replace PROBE.
CONSOLE reads positive temperature during ablation.	Thermocouple Connectors are plugged in reversed (red-to-black).	Plug Thermocouple Connectors into the matching colored CONSOLE Ports.
CONSOLE displays fault code, error code, maintenance needed, or low cylinder pressure light.	See CONSOLE User’s Manual.	

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AtriCure Inc.